

# Connecting

from

heart

to

heart



**SonoTrax**

Ultrasonic Pocket Doppler

|                         |  |
|-------------------------|--|
| Display                 | 45 x 25 mm LCD                             |
| Ratio cardiaco          | 50 ~ 210 bpm                               |
| Resolución              | 1 bpm                                      |
| Frecuencia amplificador | 4 KHz / 8 KHz                              |
| Rango grabación         | 480 s / 240 s                              |
| Audio Output            | 0,5 W                                      |
| Autodesconexión         | 1 minuto                                   |
| Medidas                 | 32 (Profundo) x 85 (Ancho) x 138 (Alto) mm |
| Peso                    | 290 ± 5 gr (bateria inc)                   |



